CTOR

Leading in Orthodontic Implants
psm-medical.com

Quality is never an accident.

It is always the result of high intention, sincere effort, intelligent direction and skillful execution; it represents the wise choice of many alternatives.

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
Introduction | Einführung

Skeletal anchorage in orthodontics.

CTOR and **BENEFIT®**. The use of **CTOR plates** in combination with our widely used palatal TAD system **BENEFIT®** considerably expands the range of orthodontic corrections. This is particularly advantageous in situations where conventional TADs are not ideally suited due to insufficient stability of the cortical bone and soft tissue or where their accessibility and practicability are limited. By utilising improved bone and soft tissue conditions without compromising the force system design, **CTOR plates** offer increased stability while providing multiple options for force and moment application.

CTOR plates are easily applied, replaced and removed without flap surgery, allowing the orthodontist to adjust the plate according to the patients' needs during the course of treatment. These adjustments can be made without surgery or anaesthesia, making these treatment options cost-effective and attractive to the patient.

The versatility and mobility of the **CTOR plates** enable the treatment of a variety of complex malocclusions without surgery combining it with the most successful **BENEFIT®** palatal TAD System.

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Developed in cooperation with:

Mani Alikhani, D.M.D, M.S., Ph.D. Harvard | School of Dental Medicine „The Forsyth Institute“

Director:

CTOR, Consortium for Transnational | Orthodontic Research | Orthodontic Research Group IADR
International Association of Dental Research

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Mani Alikhani, D. M. D. , M. S. , Ph. D. Harvard | School of Dental Medicine „Das Forsyth-Institut“

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CTOR Plates | CTOR-Platte

BENEFIT® TADs and **CTOR plates** fulfil mechanical requirements but take advantage of better cortical bone that is further away from the teeth.

BENEFIT® TADs are easy to insert and remove by the orthodontist just through the mucosa, and do not require repeated administration of anesthetics. The **CTOR plates** are connected with small fixation screws on top of the **BENEFIT**® screws.

They offer a convenient and cost-effective solution for the treatment of difficult malocclusions.

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

Features

- Insertion/removal without the need for surgical intervention
- No soft tissue contact of the **CTOR plates** – therefore no irritation
- Maximum contact with cortical bone
- Can be adjusted outside the mouth
- Easy change of force application point
- Can be used for multiple aims (dental or skeletal) simultaneously
- Multiple use of the **BENEFIT**® screws with different appliances
- Resistant to forces that loosen or destabilise conventional TADs
- Easy, on-demand replacement without additional anaesthesia
→ easy access in case of emergency

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

BENEFIT® Orthodontic Screws, Titanium

BENEFIT® KFO-Schrauben, Titan

NEW premium implants

L	Ø 2.0 mm
7 mm	ST-33-54207-2
9 mm	ST-33-54209-2
11 mm	ST-33-54211-2

P2 ST

BENEFIT® Screw Dispenser
BENEFIT® Schraubenspender

L	Ø 2.0 mm
7 mm	ST-33-54207-10
8 mm	ST-33-54208-10
9 mm	ST-33-54209-10
10 mm	ST-33-54210-10
11 mm	ST-33-54211-10

P10 ST

BENEFIT® Screws System | BENEFIT® Schrauben System

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
CTOR Plates | CTOR-Platten

DISTALIZATION | DISTALISIERUNG

Designed for posterior retraction (A). By modifying the Eagle Plate shape, and removing one wing, it can be used for unilateral distalization (B).

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


P1

RETRACTION – PROTRACTION – INTRUSION | RETRAKTION – PROTRAKTION – INTRUSION

2 T Plate types: For retraction, protraction or intrusion of posterior and anterior segments, flexible in placement in both jaws.

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


CTOR Plates | CTOR-Platten

MEZIALIZATION | MESIALISIERUNG

Wing Plates offer targeted application of force for mesialisation of maxillary teeth. Useful for protrusion of posterior teeth when anterior teeth do not require re-traction.

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
33-54780 29 mm P1


INTRUSION | INTRUSION

This Wing Plates can be used for both palatal intrusion or mesialization of dental units as shown in these schematics.

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


Instrumente & *BENEFIT*® Starter Sets

Drills | Bohrer

The instruments can be used for several **PSM Medical** products. In order to find the right instruments for your use, they are indicated in color.

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Sterile packaging | Sterile Verpackung

Implants – in single sterile packaging

The outer and inner packaging contains labels indicating UDI code-like article number, LOT number, expiration date, etc. for the patient folder and for full traceability.

The sterile screws are contained in a small colour-coded screw carrier which protects the screw and at the same time acts as a pick up tray. This carrier is placed in a clear, sealed blister package. The screw can easily be removed from the carrier by means of the self retaining screwdriver.

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


Handles | Handgriffe

The manually turned unit is connected to your existing contra-angle handpiece (blue 1:1) that enables its use like an angled screwdriver.

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Manually turned unit mod. to Pauls, with adjustable torque from 0 – 40 Ncm

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Thumb screw for Dental mandrel with limited torque (10 Ncm)

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Laboratory Screwdriver, fixed
Labor Schraubendreher, fest

33-18288		
----------	--	--

MOPS BONE Rapper | MOPS BONE Rapper

MOPS Micro Osteo Perforations

Minimally invasive acceleration of tooth movement using micro-osteoperforations (MOPs) and the BONE Rapper

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Micro Osteoperforation performed on a patient undergoing orthodontic tooth movement has been shown in a number of studies to accelerate orthodontic treatment. The micro perforations stimulate the cytokine response in the patient's alveolar bone and leads to an increased rate of bone remodeling.

The **MOPS BONE Rapper** is a sterile, minimally-invasive device that can be used in minutes by an orthodontist. It is ideal for use in conjunction with TAD's and fixed or removable treatment systems. By connecting the single-use instrument to our screwdriver handle or to a contra-angle handpiece, it is inserted under local anaesthesia into the alveolar bone in the desired area right through the soft tissue. An optional depth stop for the contra- angled version reduces the insertion depth to 5 or 7 mm. However, if desired, it can be inserted to a maximum depth of 9 mm.

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AJODO, Vol. 144, Nr. 5, Nov. 2013: Effect of Micro-Osteoperforation on the Rate of Tooth Movement

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Battery Screwdriver iSD900 | Akku Schraubendreher iSD900

Cordless Battery Screwdriver with Torque Calibration System for the contolled fixation of PSM Miniscrews.

The iSD 900 will guarantee precise torque control due to the unique – easy to use – torque calibration system (10 – 40 Ncm). Time-saving and precise, riskfree insertion and tightening of screws and abutments. The complete set contains a reliable and user friendly handpiece for Dental attachments. Reverse mode is indicated with sound. Standard rechargeable AAA Ni-Mh batteries can be used.

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BENEFIT® Starter Sets | BENEFIT® Starter Sets

BENEFIT® STARTER SET

S


33-55000	BENEFIT® STARTER SET » S	without implants ohne Implantate
-----------------	---------------------------------	------------------------------------

consisting of: | bestehend aus:

	Storage	Aufbewahrung
73-31990	BENEFIT® small instruments tray, empty	BENEFIT® kleiner Instrumentenbehälter, leer
	Instruments	Instrumente
10-67513	Drill, 1.4 x 33 mm, WL 15 mm, DENTAL, red	Bohrer, 1,4 x 33 mm, AT 15 mm, DENTAL, rot
33-18232	Screwholder for implants	Schraubenhalter für Implantate
33-18236	Thumb screw for Dental mandrel	Handgriff für Instrumente mit Dentalansatz
33-54512	Screwdriver for abutment, 20 mm	Eindrehinstrument für Abutments, 20 mm

BENEFIT® STARTER SET

M

33-55001	BENEFIT® STARTER SET » M	without implants ohne Implantate
-----------------	---------------------------------	------------------------------------

consisting of: | bestehend aus:

	Storage	Aufbewahrung
73-31990	BENEFIT® small instruments tray, empty	BENEFIT® kleiner Instrumentenbehälter, leer
	Instruments	Instrumente
10-67513	Drill, 1.4 x 33 mm, WL 15 mm, DENTAL, red	Bohrer, 1,4 x 33 mm, AT 15 mm, DENTAL, rot
33-18232	Screwholder for implants	Schraubenhalter für Implantate
33-18236	Thumb screw for Dental mandrel	Handgriff für Instrumente mit Dentalansatz
33-54512	Screwdriver for abutment, 20 mm	Eindrehinstrument für Abutments, 20 mm
33-18244	Drill Extension, + 16 mm (max. 600 rpm)	Bohrerverlängerung, + 16 mm (max. 600 U/min)
33-825-33	BENEFIT® DIRECT TX Screw holder, 28 mm	BENEFIT® DIRECT TX Schraubenhalter, 28 mm

BENEFIT® Starter Sets | BENEFIT® Starter Sets

BENEFIT® STARTER SET » L

33-55002	BENEFIT® STARTER SET » L	without implants ohne Implantate
-----------------	---------------------------------	------------------------------------

consisting of: | bestehend aus:

	Storage	Aufbewahrung
73-31990	BENEFIT® small instruments tray, empty	BENEFIT® kleiner Instrumentenbehälter, leer
	Instruments	Instrumente
10-67513	Drill, 1.4 x 33 mm, WL 15 mm, DENTAL, red	Bohrer, 1,4 x 33 mm, AT 15 mm, DENTAL, rot
33-18232	Screw holder for implants	Schraubenhalter für Implantate
33-18236	Thumb screw for Dental mandrel	Handgriff für Instrumente mit Dentalansatz
33-54512	Screwdriver for abutment, 20 mm	Eindrehinstrument für Abutments, 20 mm
33-825-33	BENEFIT® DIRECT TX Screw holder, 28 mm	BENEFIT® DIRECT TX Schraubenhalter, 28 mm
33-10903	Screw holder with stop for insertion guide	Schraubenhalter mit Stop für Insertions guide
33-54515	Screwdriver for small screws, 30 mm	Schraubendreher für kleine Schrauben, 30 mm

BENEFIT® STARTER SET » XL

33-55005	BENEFIT® STARTER SET » XL	without implants ohne Implantate
-----------------	----------------------------------	------------------------------------

consisting of: | bestehend aus:

	Storage	Aufbewahrung
73-31990	BENEFIT® small instruments tray, empty	BENEFIT® kleiner Instrumentenbehälter, leer
	Instruments	Instrumente
10-67513	Drill, 1.4 x 33 mm, WL 15 mm, DENTAL, red	Bohrer, 1,4 x 33 mm, AT 15 mm, DENTAL, rot
33-10903	Screw holder with stop for insertion guide	Schraubenhalter mit Stop für Insertions guide
33-18232	Screw holder for implants	Schraubenhalter für Implantate
33-18236	Thumb screw for Dental mandrel	Handgriff für Instrumente mit Dentalansatz
33-54512	Screwdriver for abutment, 20 mm	Eindrehinstrument für Abutments, 20 mm
33-825-33	BENEFIT® DIRECT TX Screw holder, 28 mm	BENEFIT® DIRECT TX Schraubenhalter, 28 mm
33-18244	Drill Extension, + 16 mm (max. 600rpm)	Bohrerverlängerung, + 16 mm (max. 600 U/min)
33-18288	Laboratory screwdriver, fixed	Labor Schraubendreher, fixiert

Catalogues | Kataloge

More innovative **PSM** products are presented in our other catalogues. Please request it from your local **PSM** distribution partner or download from our website. www.psm-medical.com.

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BENEFIT

With the introduction of the **BENEFIT**® System **PSM** has considerably enhanced the therapeutic diversity of the treatments available for Mesialization, Distalization, RPE as well as uprighting or anchorage of the dentition. With **BENEFIT**®, **PSM** was the first company to offer a full range of skeletal anchorage appliances for the palate – a trend that has become „state of the art“ in Orthodontics.

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smartlock

The new **PSM** – **smartlock**® system is largely based on the successful and proven **BENEFIT**® system. **smartlock**® focuses on the essentials: Unique comfort as well as safe and easy handling that doctors and patients can rely on. The trend-setting concept, combined with the best materials and the most modern manufacturing processes, guarantees a treatment result at the highest level.

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quattro

The **quattro**® Mini Anchor System is known as one of the first and most successful orthodontic Mini-Implant Systems worldwide. Based on this system the **BENEFIT**® screws were developed. This enables the user to utilize the entire scope of orthodontic implantation with unmatched possibilities and compatible instrumentation.

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mentoplate

The **mentoplate**® system provides a perfect solution for enhanced stability with special orthodontic bone anchored plates for class III correction. The System includes Titanium plates which are fixed to the bony surface with selfdrilling screws.

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TADMAN®

The complete digital workflow realized with no efforts for the orthodontist. Based on the **BENEFIT**® System. For more information: www.tadman.de

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Specifications and design are subject to change without prior notice. All pictures, photos and product descriptions are the property of PSM Medical GmbH. Utilization and copies by third parties have to be authorized by PSM Medical GmbH. All rights reserved! All trademarks and brands mentioned in this brochure are subject to unrestricted regulations of the trade mark laws valid and right of ownership of the registered owner.

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PSM North America, Inc.